## To find a Dealer near you contact the MEDA Office @ 608-588-7878 / medaoffice@chemstarworks.com





The Power of Chlorine Dioxide in an economical, convenient concentrate...

#### PREMIUM PROTECTON

- CHLORINE DIOXIDE Primary germicide that kills all major mastitis causing organisms FAST including Staph, Strep, Mycoplasma, molds & yeasts.
- LAS -Secondary Germicide that provides germicidal persistency ensuring teats are protected right up to the next milking.
- EFFECIVE TEAT MARKING Verify teats are being properly dipped to milker to milker and milking to milking.

#### Express for Pre-Dipping

- Cleans & Kills quickly & efficiently
- 15 second Kill
- Yellow or Blue Marking for easy ID
- Residue Safe

### Express for Post Dipping

- Fast, Broad Spectrum Complete Kill
  - · Chlorine Dioxide Germicide
- Extends Kill / Protection until next Milking
  - LAS Germicide
- Yellow or Blue Marking of easy ID
- Emollients = Propylene Glycol & Lanolin

#### AOAC Sanitizing Efficacy Results—

6 log, 15 Second Kill under 10% Milk Load.
(A Log Reduction of 6 represents 99.9999% Reduction in Mastitis
Causing Organisms)

# DUAL Germicidal ACTION!





**DIP - FOAM - SPRAY** 



Ask us about...
DISPENSING OPTIONS!

Message:





# **Express**®

The Power of Chlorine Dioxide in an economical, convenient concentrate...

Ask us about...

DISPENSING OPTIONS!



The AccuMAX accurately dilutes dip at the push of a button. No electricity is needed... just a water outlet (20 psi). Water moves through the venturi pump system pulling the correct amount of concentrated product through metering tips to make RTU Product. It is great for filling dip cups or gallon jugs.



The Parlor Master accurately dilutes and delivers product to the parlor with one convenient unit. No electricity is needed... just a water outlet (20 psi) and compressed air. This pump has a float inside and regenerates product as it is depleted, leaving fresh product. Metering tips are used to insure the right amount of concentrated product is used.



The Knight system has a micro processor based digital blending controller which allows you to use different teat dip concentrations. It has high reliable electric diaphragm pumps which allow it to deliver precise fluid volumes based on flow meter signals and formulation composition. It also has a SD card which allows you to retrieve usage reports.

### Message: